

Human CD38 protein (Active Recombinant) (C-His) (STJP000926)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	952
Gene Symbol	CD38
UniProt ID	CD38_HUMAN
Immunogen Region	Val43-Ile300
Immunogen Sequence	VPRWRQQWSPGTTKRFPEV VLARCVKYTEIHPEMRHVDC QSVWDAFKGAFISKHPCNIT EEDYQPLMKLGTTQVPCNKI LLWSRIKDLAQFTQVQRDM FTLEDTLGLYADDLTWCGE FNTSKINYQSCPDWRKDCSN NPVSFVWKTVSRFAEAACD VVHVMLNGSRSKIFDKNSTF GSVEVHNLQPEKVQTLAWV IHGGREDSRDLCDPTIKEL ESIISKRNQFSCCKNIYRP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

