

Human CD5 protein (Recombinant) (C-His) (STJP000933)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
UniProt ID	CD5_HUMAN
Immunogen Region	Arg25-Asn371
Immunogen Sequence	RLSWYDPDFQARLTRNSKQ QGQLEVYLKDGWHMVCSSQSW GRSSKQWEDPSQASKVCQRL NCGVPLSLGPFVLYTPQSS IICYGQLGSFNSCSHSRNDM CHSLGLTCLPQKTTPTTR PPPTTTPPTAPPRLQLVAQ SGGQHCAGVVEFYSGSLGGT ISYEAQDKTQDLENFLCNNL QCGSFLKHLPETEAGRAQDP GEPREHQPLPIQWKIQNSSC TSLEHCFRKKIKPKSGRVL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

