

Human Siglec-3/CD33 protein (Recombinant) (C-hFc&His) (STJP000934)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	945
Gene Symbol	CD33
UniProt ID	CD33_HUMAN
Immunogen Region	Asp18-His259
Immunogen Sequence	DPNFWLQVQESVTVQEGLCV LVPCTFFHPIPYDKNPVS GYWFREGAIIISRDSPVATNK LDQEVQEETQGRFRLLDPS RNNCSLSIVDARRRDNGSYF FRMERGSTKYKSPQLSVH VTDLTHRPKILIPGTLEPGH SKNLTCSVSWACEQGTPPIF SWLSAAPTSLGPRTHSSVL IITPRPQDHGTNLTQVKFA GAGVTERTIQLNVTYVQPN PTTGIFPGDGGSGKQETRAG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

