

Human C1qR/CD93 protein (Recombinant) (C-hFc&His) (STJP000939)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22918
Gene Symbol	CD93
UniProt ID	C1QR1_HUMAN
Immunogen Region	Ala24-Lys580
Immunogen Sequence	ADTEAVVCVGTACYTAHSGK LSAAEAQNHCNQNGGNLATV KSKEEAQHVQRVLAQLLRRE AALTARMSKFWIGLQREK GK CLDPSLPLKGF SWVGGG EDT PYSNWHKELRNSCISKRCVS LLLDLSQPLLPSRLPKWSEG PCGSPGSPGSNIEGFVCKFS FKGMCRLALGGPGQVTTYT PFQTTSSSLEAVPFASAANV ACGEGDKDETQSHYFLCKEK APDVF DWGSSGPLCVSPKY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

