

Human JAM-3 (A149P) protein (Recombinant) (C-hFc&His) (STJP000945)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83700
Gene Symbol	JAM3
UniProt ID	JAM3_HUMAN
Immunogen Region	Val32-Asn241 Ala149Pro
Immunogen Sequence	VNLKSSNRTPVVQEFESVEL SCIITDSQTSDPRIEWKKIQ DEQTTYVFFDNKIQGDLAGR AEILGKTSLSKIWNVTRRDSA LYRCEVVARNDRKEIDEIVI ELTVQVKPVTVPVCRVPKVPV VGKMATLHCQESEGHPRPY SWYRNDVPLPTDSRANPRFR NSSFHLNSETGLVFTAVHK DDSGQYYCIASNDAGSARCE EQEMEVDLN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

