

Human MIC-A protein (Recombinant) (C-His) (STJP000947)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100507436
Gene Symbol	MICA
UniProt ID	MICA_HUMAN
Immunogen Region	Met1-Gln308
Immunogen Sequence	MGLGPVFLLAGIFPFAPPG AAAPHSRLRYNLTVLSWDGS VQSGFLTEVHLDGQPFLRCD RQKCRAPQGQWAEDVLGNK TWDRETRDLTGNGKDLRMTL AHKDKQKEGLHSLQEIRVCE IHEDNSTRSSQHFYYDGELF LSONLETKEWTMPQSSRAQT LAMNVRNFLKEDAMKTKTHY HAMHADCLQLRRYLKSGVV LRRTVPPMVNVTREASEGN ITVTCRASGFYPWNITLSW

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

