

Human IL-12A protein (Recombinant) (C-hFc&His) (STJP000948)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3592
Gene Symbol	IL12A
UniProt ID	IL12A_HUMAN
Immunogen Region	Arg23-Ser219
Immunogen Sequence	RNLPVATPDPGMFPCPLHHSQ NLLRAVSNMLQKARQTLEFY PCTSEEIDHEDITKDKTSTV EACLPLELTKNESCLNSRET SFITNGSCLASRKTSFMMAL CLSSIYEDLKMYQVEFKTMN AKLLMDPKRQIFLDQNMLAV IDELMQALNFNSETVPQKSS LEEPDFYKTIKILLCILLHAF RIRAVTIDRVMSYLNAS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

