

Human Cathepsin S protein (Active Recombinant) (C-His) (STJP000955)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50mM MES, 100mM NaCl, pH 6.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1520
Gene Symbol	CTSS
UniProt ID	CATS_HUMAN
Immunogen Region	Gln17-Ile331
Immunogen Sequence	QLHKDPTLDHHWHLWKITYG KQYKEKNEEAVRRLIWEKNL KFMLHNLHSMGMHSYDLG MNHLDGMTSEEVMSSLR VPSQWQRNITYKSNPNRILP DSVDWREKGCVTEVKYQGSC GACWAFSAVGALEAQLKLT GKLVSLSAQNLDVDCSTEKYG NKGCGGGFMTTAFQYIIDNK GIDSDASYPYKAMDQKCQYD SKYRAATCSKYTELPYGRD VLKEAVANKGPVSVGVDR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

