

Human IGFBP-4 protein (Recombinant) (C-His) (STJP000959)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3487
Gene Symbol	IGFBP4
UniProt ID	IBP4_HUMAN
Immunogen Region	Met1-Glu258
Immunogen Sequence	MLPLCLVAALLAAGPGPSL GDEAIHCPPCSEELARCRP PVGCEELVREPGCGCCATCA LGLGMPGCVYTPRCGSLRC YPPRGVEKPLHTLMHGQGVG MELAEIEAIQESLQPSDKDE GDHPNNSFSPCSAHDRRCLQ KHFAKIRDRTSGGKMKVNG APREDARPVPQGSCQSELHR ALERLAASQSRTHEDLYIIP IPNCDRNGNFHPKQCHPALD GQRGKCWCVDRTGKVLPG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

