

Human WIF-1 (Q166K) protein (Recombinant) (C-His) (STJP000963)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM HAc-NaAc, 150mM NaCl, 0.5% CHAPS, pH 6.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11197
Gene Symbol	WIF1
UniProt ID	WIF1_HUMAN
Immunogen Region	Gly29-Trp379Gln166Lys
Immunogen Sequence	GPPQEESLYLWIDAHQARVL IGFEEDILIVSEGKMAPFTH DFRKAQQRMPAIPVNIHSMN FTWQAAGQAEYFYEFSLRS LDKGIMADPTVNVPLLGTVP HKASVVQVGFPCLGKQDGV AFEVDVIVMNSEGNTILKTP QNAIFFKTCQQAECPGGCRN GGFCNERRICECPDGFHGPH CEKALCTPRCMNGGLCVTPG FCICPPGFYGVNCDKANCST TCFNGGTCTFYPGKCICPPG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

