

Human Contactin 2/CNTN2 protein (Recombinant) (C-His) (STJP000964)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6900
Gene Symbol	CNTN2
UniProt ID	CNTN2_HUMAN
Immunogen Region	Ser31-Asn1012
Immunogen Sequence	SQTTFGPVPFEDQPLSVLFPE ESTEEQVLLACRARASPPAT YRWKMNGTEMKLEPGSRHQL VGGNLVIMNPTKAQDAGVYQ CLASNPVGTVVSREAILRFG FLQEFKSKEERDPVKAHEGWG VMLPCNPPAHYPGLSYRWLL NEFPNFIPTDGRHFVSQTTG NLYIARTNASDLGNYSLAT SHMDFSTKSVFSKFAQLNLA AEDTRLFAPSIKARFPAETY ALVGQQVTLECFAGNPVP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

