

Human Fetuin A/AHSG protein (Recombinant) (C-His) (STJP000968)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	197
Gene Symbol	AHSG
UniProt ID	FETUA_HUMAN
Immunogen Region	Ala19-Val367
Immunogen Sequence	APHGPGLIYRQPNCDDEPETE EALVAIDYINQNLPGWGYKH TLNQIDEVKVWPQQPSGELF EIEIDTLETTCHVLDPTPVA RCSVRQLKEHAVEGDCDFQL LKLDGKFSSVYAKCDSSPDS AEDVRKVCQDCPLLAPLNDT RRVHAAKAALAAFNQNNNGS NFQLEEISRAQLVPLPPSTY VEFTVSGTDCVAKEATEAAK CNLLAEKQYGFKATLSEKL GGAEVAVTCMVFTQTPVSS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

