

Human GDNFR-alpha-2/GFRA2 protein (Recombinant) (C-His) (STJP000972)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2675
Gene Symbol	GFRA2
UniProt ID	GFRA2_HUMAN
Immunogen Region	Ser22-Ser441
Immunogen Sequence	SPSSLQGPELHGWRPPVDCV RANELCAAESNCSSRYRTL R QCLAGRDRNTMLANKECQAA LEVLQESPLYDCRCKRGMKK ELQCLQIYWSIHLGLTEGEE FYEASPYEPVTSRLSDIFRL ASIFSGTGADPVVSAKSNHC LDAAKACNLNDNCKKLRSY ISICNREISPTERCNRRKCH KALRQFFDRVPSEYTYRMLF CSCQDQACAERRRQTILPSC SYEDKEKPNCLDLRGVCRT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

