

Human IL-6 protein (Active Recombinant) (C-His) (STJP000974)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
UniProt ID	IL6_HUMAN
Immunogen Region	Val30-Met212
Immunogen Sequence	VPPGEDSKVAAPHRQPLTS SERIDKQIRYILDGIALRK ETCNKSNMCESSKEALAENN LNLPKMAEKDGCFSQGFNEE TCLVKIITGLLEFEVYLEYL QNRFESSEEQARAVQMSTKV LIQFLQKKAKNLDAITTPDP TTNASLLTKLQAQNLWLQDM TTHLILRSFKEFLQSSLRAL RQM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

