

Human Frizzled-7/FZD7 protein (Recombinant) (C-hFc&His) (STJP000977)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8324
Gene Symbol	FZD7
UniProt ID	FZD7_HUMAN
Immunogen Region	Gln33-Leu185
Immunogen Sequence	QPYHGEKGISVPDHGFCQPI SIPLCTDIAYNQITLPNLLG HTNQEDAGLEVHQFYPLVKV QCSPELRFFLCSMYAPVCTV LDQAIPPCRSLCERARQGCE ALMNKFGFQWPERLRCENFP VHGAGEICVQGNTSDGSGGP GGGPTAYPTAPYL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

