

Human CD301 protein (Recombinant) (C-His) (STJP000979)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10462
Gene Symbol	CLEC10A
UniProt ID	CLC10_HUMAN
Immunogen Region	Gln61-His316
Immunogen Sequence	QNSKFQRDLVTLRTDFSNT SNTVAEIQALTSQGSSEET IASLKAEEVEGFKQERQAGVS ELQEHTTQKAHLGHCPHCPS VCVPVHSEMLLRVQQLVQDL KKLTCQVATLNNNASTEGTC CPVNWVEHQDSCYWFSHSGM SWAEAEKYCQLKNAHLVVIN SREEQN FVQKYLGSAYTWMG LSDPEGAWKWVDGTDYATGF QNWKPGQPDWQGHGLGGGE DCAHFHPDGRWNDDVCQRP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

