

Human Dkk-3 protein (Recombinant) (C-His) (STJP000981)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27122
Gene Symbol	DKK3
UniProt ID	DKK3_HUMAN
Immunogen Region	Ala22-Ile350
Immunogen Sequence	APAPTATSAPVKPGPALSY QEEATLNEMFREVEELMEDT QHKLRSVEEMEAEAAAAKA SSEVNLANLPSPYHNETNTD TKVGNNTIHVHREIHKITNN QTGQMVFSETVITSVGDEEG RRSHECIIDEDCGPSMYCQF ASFQYTCQPCRGQRMLCTRD SECCGDQLCVWGHCTKMATR GSNGTICDNQRDCQPLCCA FQRGLLPVCTPLPVEGELC HDPASRLDLITWELEPDG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

