

Human MMP-1 protein (Recombinant) (C-His) (STJP000983)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4312
Gene Symbol	MMP1
UniProt ID	MMP1_HUMAN
Immunogen Region	Phe20-Asn469
Immunogen Sequence	FPATLETQEQDVLVQKYLE KYYNLKNDRQVEKRRNSGP VVEKLGKMQEFFFGLKVTGKP DAETLKVMKQPRCGVPDVAQ FVLTEGNPRWEQTHLYRIE NYTPDLPRADVDHAIKAFQ LWSNVTPLTFTKVSEGGADI MISFVRGDHRDNSPFDGPGG NLAHAFQPGPGIGGDAHFDE DERWTNNFREYNLHRVAAHE LGHSLGLSHSTDIGALMYPYS YTFSGDVQLAQDDIDGIQA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

