

Human Cathepsin D protein (Recombinant) (C-His) (STJP000985)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1509
Gene Symbol	CTSD
UniProt ID	CATD_HUMAN
Immunogen Region	Leu21-Leu412
Immunogen Sequence	LVRIPLHKFTSIRRTMSEVG GSVEDLIAKGPVSKYSQAVP AVTEGPIPEVLKNYMDAQYY GEIGIGTPPQCFTVVFDTS SNLWVPSIHCKLLDIACWIH HKYNSDKSSTYVKNGTSFDI HYGSGSLSGYLSQDTVSVPC QSASSASALGGVKVERQVFG EATKQPGITFIAAKFDGILG MAYPRISVNNVLPVFDNLMQ QKLVDQNIFSFLSRDPDAQ PGGELMLGGTDSKYYKGSL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

