

Human IGFBP-2 protein (Recombinant) (C-His) (STJP000986)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3485
Gene Symbol	IGFBP2
UniProt ID	IBP2_HUMAN
Immunogen Region	Glu40-Gln328
Immunogen Sequence	EVLFRCPPCTPERLAACGPP PVAPPAVA AVAGGARMPCA ELVREPGCGCCSVCARLEGE ACGVYTPRCGQGLRCYPHPG SELPLQALVMGEGTCEKRRD AEYGASPEQVADNGDDHSEG GLVENHVDSTMNMLGGGGSA GRKPLKSGMKELAVFREKVT EQHRQMGKGGKHHLGLEPK KLRPPPARTPCQQLDQVLE RISTMRLPDERGPLEHLYSL HIPNCDKHGLYNLKQCKMS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

