

Human Osteonectin/SPARC protein (Recombinant) (C-His) (STJP000988)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6678
Gene Symbol	SPARC
UniProt ID	SPRC_HUMAN
Immunogen Region	Ala18-Ile303
Immunogen Sequence	APQQEALPDETEVVEETVAE VTEVSVGANPVQVEVGEFDD GAEETEEVVAENPCQNHHC KHGKVCELDENNTPMCVCQD PTSCPAPIGEFEKVCNSDNK TFDSSCHFFATKCTLEGTKK GHKLHLDYIGPCKYIPPCLD SELTEFPLRMRDWLKNVLVT LYERDEDNNLLTEKQKLRVK KIHENEKRLEAGDHPVELLA RDFEKNYNMYIFPVHWQFGQ LDQHPIDGYLSHTELAPLR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

