

Human B7-H4/VTCN1 protein (Recombinant) (C-His) (STJP000989)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79679
Gene Symbol	VTCN1
UniProt ID	VTCN1_HUMAN
Immunogen Region	Phe29-Ala258
Immunogen Sequence	FGISGRHSITVTTVASAGNI GEDGILSCTFEPDIKLSDIV IQWLKEGVLGLVHEFKEGKD ELSEQDEMFRGRTAVFADQV IVGNASRLKKNVQLTDAGTY KCYIITSKGKGNANLEYKTG AFSMPEVNVVDYNASSETLRC EAPRWFPQPTVVWASQVDQG ANFSEVSNTSFELNSENVTM KVVSVLYNVTINNTYSCMIE NDIAKATGDIKVTSEIKRR SHLQLLSKA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

