

Human MMP-3 protein (Active Recombinant) (C-His) (STJP000991)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4314
Gene Symbol	MMP3
UniProt ID	MMP3_HUMAN
Immunogen Region	Tyr18-Cys477
Immunogen Sequence	YPLDGAARGEDTSMNLVQKY LENYYDLEKDVKQFVRRKDS GPVVKKIREMQKFLGLEVTG KLDSDTLLEVMRKPRCGVPDV GHFRTFPGIPKWRKTHLYR IVNYTPDLPKDAVDSAVEKA LKVWEEVTPLTFSRLYEGEA DIMISFAVREHGDYFPDGP GNVLAHAYAPGPGINGDAHF DDDEQWTKDTTGTNLFLVAA HEIGHSLGLFHSANTEALMY PLYHSLTDLTRFRLSQDDI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

