

Human MMP-13 protein (Recombinant) (C-His) (STJP000993)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4322
Gene Symbol	MMP13
UniProt ID	MMP13_HUMAN
Immunogen Region	Leu20-Cys471
Immunogen Sequence	LPLSGGDEDDLSEEDLQFA ERYLRSYYHPTNLAGILKEN AASSMTERLREMQSFFGLEV TGKLDNDNTLDVMKKPRCGVP DVGEYNVFPRTLKWSKMNLT YRIVNYTPDMTHSEVEKAFK KAFKVVSDVTPLNFTRLHDG IADIMISFGIKEHGDFYFPD GPSGLLAHAFFPPGPNYGGDA HFDDDETWTSSSKGYNLFLV AAHEFGHSLGLDHSKDPGAL MFPIYTYTGKSHFMLPDDD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

