

Human Kallikrein-8/KLK8 protein (Recombinant) (C-His) (STJP000997)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11202
Gene Symbol	KLK8
UniProt ID	KLK8_HUMAN
Immunogen Region	Gln29-Gly260
Immunogen Sequence	QEDKVLGGHECQPHSQPWQA ALFQGQQLLGGVVLVGGNWV LTAAHCKKPKYTVRLGDHSL QNKDGPEQEIPVVQSIPHC YNSSDVEDHNDLMLLQLRD QASLGSKVKPISLADHCTQP GQKCTVSGWGTVTSPRENFP DTLNCAEVKIFPQKKCEDAY PGQITDGMVCAGSSKGADTC QGDSGGPLVCDGALQGITSW GSDPCGRSDKPGVYTNICRY LDWIKKIIGSKG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

