

Human Hedgehog-interacting protein/HHIP protein (Recombinant) (C-His) (STJP000999)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64399
Gene Symbol	HHIP
UniProt ID	HHIP_HUMAN
Immunogen Region	Phe18-Asp670
Immunogen Sequence	FFEGDAKFGERNEGSGARRR RCLNGNPPKRLKRRDRRMMS QLELLSGGEMLCGGFYPRLS CCLRSDSPGLGRLENKIFSV TNNTECGKLLLEEIKCALCSP HSQSLFHSPEREVLERDLVL PLLCKDYCKEFFYTCTRGHIP GFLQTTADEFCFYARKDGG LCFPDFPRKQVRGPASNYLD QMEEYDKVEEISRKHKNCF CIQEVVSGLRQPVGALHSGD GSQRLFILEKEGYVKILTP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

