

Human IFN-lambda 2/IL-28A protein (Recombinant) (C-His) (STJP001000)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	282616
Gene Symbol	IFNL2
UniProt ID	IFNL2_HUMAN
Immunogen Region	Val26-Val200
Immunogen Sequence	VPVARLHGALPDARGCHIAQ FKSLSPELQAFKRAKDALE ESKLLKDCRCHSRFPRTWD LRQLQVREPRMALEAELALT LKVLEATADTDPALVDVLDQ PLHTLHHILSQFRACIQQP TAGPRTRGRLLHHLWYRLQEA PKKESPGCLEASVTFNLFRL LTRDLNCVASGDLCV

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

