

Human DNER protein (Recombinant) (C-hFc&His) (STJP001002)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92737
Gene Symbol	DNER
UniProt ID	DNER_HUMAN
Immunogen Region	Arg29-His637
Immunogen Sequence	RGSSLANPVPAAPLSAPGPC AAQPCRNGGVCTSRPEPDQ HPAPAGEPGYSCTCPAGISG ANCQLVADPCASNPCHHGNC SSSSSSSSDGYLCICNEGYE GPNCEQALPSLPATGWTESM APRQLQVPVATQEPDKILPR SQATVTLPTWQPKTGQKVVE MKWDQVEVIPDIACGNASSN SSAGGRLVSFEVPQNTSVKI RQDATASLILLWVKTATGFQ QCSLIDGRSVTPLQASGGL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

