

Human TIM-3/HAVCR2/CD366 protein (Recombinant) (C-hFc&His) (STJP001003)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84868
Gene Symbol	HAVCR2
UniProt ID	HAVR2_HUMAN
Immunogen Region	Ser22-Arg200
Immunogen Sequence	SEVEYRAEVGQNAYLPCFYT PAAPGNLVPVCWGKGACPVF ECGNVVLRTDERDVNYWTSR YWLNGDFRKGDVSLTIENVT LADSGIYCCRIQIPGIMNDE KFNLLKLVIPAKVTPAPTRQ RDFTAAFPRLTTRGHGPAE TQTLGSLPDINLTQISTLAN ELRDSRLANDLRD SGATIR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

