

## Human IGFBP-1 protein (Recombinant) (C-His) (STJP001006)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                 |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">3484</a>                                                                                                                                                                                                                                            |
| Gene Symbol        | <a href="#">IGFBP1</a>                                                                                                                                                                                                                                          |
| UniProt ID         | <a href="#">IBP1_HUMAN</a>                                                                                                                                                                                                                                      |
| Immunogen Region   | Ala26-Asn259                                                                                                                                                                                                                                                    |
| Immunogen Sequence | APWQCAPCSAEKLALCPPVS ASCSEVTRSAGCGCCPMCAL PLGAACGVATARCARGLSCR<br>ALPGEQQPLHALTRGQGACV QESDASAPHAEEAGSPESPE STEITEEELLDNFHLMAPSE<br>EDHSILWDAISTYDGSKALH VTNIKKWKEPCRIELYRVVE SLAKAQETSGEEISKFYLPN<br>CNKNGFYHSRQCETSMDGEA GLCWCVPYPWNGKRIPGSPEI RGDPCNQIYFNVQN |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

