

Human CD2 protein (Recombinant) (C-His) (STJP001010)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	914
Gene Symbol	CD2
UniProt ID	CD2_HUMAN
Immunogen Region	Lys 25-Asp 209
Immunogen Sequence	KEITNALETWGALGQDINLD IPSFQMSDDIDDIKWEKTS D KKKIAQFRKEKETFK EKEDTY KLFKNGTLKIKHLKTDQDI YKVSIDYTKGKNVLEKIFDL KIQERVSKPKISWTCINTTL TCEVMNGTDPELNLYQDGKH LKLSQRVITHKWTTSLSAKF KCTAGNKVSKESSVEPVS CP EKGLD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

