

Human CNTFR-alpha protein (Recombinant) (C-His) (STJP001011)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1271
Gene Symbol	CNTFR
UniProt ID	CNTFR_HUMAN
Immunogen Region	Gln 23-Pro 346
Immunogen Sequence	QRHSPQEAPHVQYERLGSDV TLPCGTANWDAAVTWVRVNGT DLAPDLLNGSQLVLHGLELG HSGLYACFHRDSWHLRHQVL LHVGLPPREPVLSCRSNTYP KGFYCSWHLPTPTYIPNTFN VTVLHGSKIMVCEKDPALKN RCHIRYMHLEFSTIKYKVSIS VSNALGHNATAITFDEFTIV KPDPPENVVARPVPSNPRRL EVTWQTPSTWPDPEFPLKF FLRYRPLLDQWQHVELSD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

