

Human Ephrin-B1/EFNB1 protein (Recombinant) (C-hFc&His) (STJP001018)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1947
Gene Symbol	EFNB1
UniProt ID	EFNB1_HUMAN
Immunogen Region	Leu 28-Gly 232
Immunogen Sequence	LAKNLEPVSWSLNPFLSG KGLVIYPKIGDKLDIICPRA EAGRPYEYYKLYLVRPEQAA ACSTVLDPNVLVTCNRPEQE IRFTIKFQEFSPNYMGLEFK KHHDYITSTSNGLSLEGLN REGGVCRTRTMKIIMKVGQD PNAVTPEQLTTSRPSKEADN TVKMATQAPGSRGSLGSDSG KHETVNQEEKSGPGASGGSS GDPDG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

