

Human Chitinase-1/CHIT1 protein (Recombinant) (C-His) (STJP001019)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1118
Gene Symbol	CHIT1
UniProt ID	CHIT1_HUMAN
Immunogen Region	Ala22-Asn466
Immunogen Sequence	AKLVICYFTNWAQYRQGEARF LPKDLDP LSLCTH LIYAFAGM TNHQLSTTEW NDETLYQEFN GLKKMNP K LK TLLA IGGWNF STQKFTDMVATANNRQTFVN SAIRFLRKYSFDGLDL DWEY PGSQGSPAVDKERFTTLVQD LANAFQQE AQTSGKERLLLS AAVPAGQTYVDAGYEVDKIA QNLDFVNL MAYDFHGSWEKV TGHNSPLYKRQEESGAAASL NVDAAVQQWLQK GTPASKL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

