

Human Arginase 1/ARG1 protein (Recombinant) (C-His) (STJP001020)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	383
Gene Symbol	ARG1
UniProt ID	ARG1_HUMAN
Immunogen Region	Met 1-Lys 322
Immunogen Sequence	MSAKSRTIGIIGAPFSKGQP RGGVEEGPTVLRKAGLLEKL KEQECDVKDYGDLPFADIPN DSPFQIVKNPRSVGKASEQL AGKVAEVKKNGRISLVGGD HSLAIGSISGHARVHPDLGV IWVDAHTDINTPLTTTSGNL HGQPVSFLLKELKGKIPDVP GFSWVTPCISAKDIVYIGLR DVDPGEHYILKTLGIKYFSM TEVDRLGIGKVMEEETSYLL GRKKRPIHLSFDVDGLDPS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

