

Human Serum albumin/ALB protein (Recombinant) (C-His) (STJP001024)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	213
Gene Symbol	ALB
UniProt ID	ALBU_HUMAN
Immunogen Region	Asp 25-Leu 609
Immunogen Sequence	DAHKSEVAHRFKDLGEENFK ALVLIAFAQYLQQCPFEDHV KLVNEVTEFAKTCVADESAE NCDKSLHTLFGDKLCTVATL RETYGEMADCCAKQEPEPNE CFLQHKDDNPNLPRLVPRPEV DVMCTAFHDNEETFLKKYLY EIARRHPYFYAPELLFFAKR YKAAFTECCQAADKAACLLP KLDEL RDEGKASSAKQRLKC ASLQKFGERAFAKAWAVARLS QRFPKAFAEVSKLVTDLT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

