

## Human LILRB2/ILT-4/CD85d protein (Recombinant) (C-His) (STJP001032)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                     |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">10288</a>                                                                                                                                                                                                                                               |
| Gene Symbol        | <a href="#">LILRB2</a>                                                                                                                                                                                                                                              |
| UniProt ID         | <a href="#">LILRB2_HUMAN</a>                                                                                                                                                                                                                                        |
| Immunogen Region   | Gln 22-Val 461                                                                                                                                                                                                                                                      |
| Immunogen Sequence | QTGTIPKPTLWAEPSVITQ GSPVTLSCQGSLEAQEYRLY REKKSASWITRIRPELVKNG<br>QFHIPSITWEHTGRYGCQYY SRARWSELSDPLVLVMTGAY PKPTLSAQPSPVVTSGGRVT<br>LQCESQVAFGGFILCKEGED EHPQCLNSQPHARGSSRAIF SVGPVSPNRRWSHRCYGYDL<br>NSPYVWSSPSDLLELLVPGV SKKPSLSVQPGPVVAPGESL TLQCVSDVG YDRFVLYKEG |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

