

Human ICAM-1/CD54 protein (Recombinant) (C-His) (STJP001041)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3383
Gene Symbol	ICAM1
UniProt ID	ICAM1_HUMAN
Immunogen Region	Gln28-Glu480
Immunogen Sequence	QTSVSPSKVILPRGGSVLVT CSTSCDQPKLLGIETPLPKK ELLLPGNNRKVYELSNVQED SQPMCYSNCPDQGQSTAKTFL TVYWTPERVELAPLPSWQPV GKNLTLRCQVEGGAPRANLT VLLRGEKELKREPAVGPEA EVTTTTLVRRDHHGANFSCR TELDLRPQGLELFENTSAPY QLQTFVLPATPPQLVSPRVL EVDTQGTVVCSLDGLFPVSE AQVHLALGDQRLNPTVTYG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

