

Human NKAT-2/KIR2DL3/CD158b2 protein (Recombinant) (C-hFc&His) (STJP001044)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3804
Gene Symbol	KIR2DL3
UniProt ID	K12L3_HUMAN
Immunogen Region	His22-His245
Immunogen Sequence	HEGVHRKPSLLAHPGPLVKS EETVILQCWSDVRFQHFLLH REGKFKDTLHLIGEHHDGVS KANFSIGPMMDLAGTYRCY GSVTHSPYQLSAPSDPLDIV ITGLYEKPSLSAQPGPTVLA GESVTLSCSSRSSYDMYHLS REGEAHERRFSAGPKVNGTF QADFPLGPATHGGTYRCFGS FRDSPYEWSNSSDPLLVSVT GNPSNSWLSPTEPSSETGNP RHLH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

