

Human JAM-A/F11R/CD321 protein (Recombinant) (C-His) (STJP001046)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
UniProt ID	JAM1_HUMAN
Immunogen Region	Ser28-Ala242
Immunogen Sequence	SVTVHSSEPEVRIPENNPVK LSCAYSGFSSPRVEWKFDQG DTTRLVCYNNKITASYEDRV TFLPTGITFKSVTREDTGT TCMVSEEGGNSYGEVKVKLI VLVPPSKPTVNPSSATIGN RAVLTCSQDGSPPSEYTW F KDGIVMPTNPKSTRAFSNSS YVLNPTTGELVFDPLSASDT GEYSCEARNGYGTPMSTNAV RMEAVERNVGVIVAA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

