

Human Carboxypeptidase E/CPE protein (Recombinant) (C-His) (STJP001051)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1363
Gene Symbol	CPE
UniProt ID	CBPE_HUMAN
Immunogen Region	Met 1-Ser 453
Immunogen Sequence	MAGRGG SALLALCGALAACG WLLGAEAEQEPGAPAAGMRRR RRLQQEDGISFEYHRYPELR EALVSVWLQCTAISRIYTVG RSFEGRELLVIELSDNPGVH EPGEPEFKYIGNMHGNEAVG RELLIFLAQYLCNEYQKGNE TIVNLIHSTRIHIMP SLNPD GFEKAASQPGELKDWVFGRS NAQGIDLNRNFPDLDRIVYV NEKEGGPNNHLLKNMKKIVD QNTKLAPETKAVIHWIMDI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

