

Human uPAR/PLAUR/CD87 protein (Recombinant) (C-His) (STJP001063)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5329
Gene Symbol	PLAUR
UniProt ID	UPAR_HUMAN
Immunogen Region	Leu23-Arg303
Immunogen Sequence	LRCMQCKTNGDCRVEECALG QDLCRTTIVRLWEEGEELEL VEKSCETHSEKTNRTLSTYRTG LKITSLTEVVCGLDLCNQGN SGRAVTYSRSRYLECI SCGS SDMSCERGRHQLQCRSPPEE QCLDVVTHWQEGEEGRPKD DRHLRGCGYLPGCPGNGFH NNDTFHFLKCCNTTKCNEGP ILELENLPQNGRQCYSCKGN STHGCSSEETFLIDCRGPMN QCLVATGTHEPKNQSYMVR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

