

Human Transthyretin/TTR protein (Recombinant) (C-His) (STJP001067)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7276
Gene Symbol	TTR
UniProt ID	TTHY_HUMAN
Immunogen Region	Gly 21-Glu 147
Immunogen Sequence	GPTGTGESKCPLMVKLVDAV RGSPAINVAVHVFRKAADD T WEPFASGKTSESGELHGLTT EEEFVEGIYKVEIDTKSYWK ALGISPFHEHAEVVFTANDS GPRRYTIAALLSPYSYSTTA VVTNPKE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

