

Human JAM-A/F11R/CD321 protein (Recombinant) (C-hFc&His) (STJP001068)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
UniProt ID	JAM1_HUMAN
Immunogen Region	Ser28-Ala242
Immunogen Sequence	SVTVHSSEPEVRIPENNPVK LSCAYSGFSSPRVEWKFDDG DTTRLVCYNNKITASYEDRV TFLPTGITFKSVTREDTGTY TCMVSEEGGNSYGEVKVKLI VLVPPSKPTVNPSSATIGN RAVLTCSAQDGGSPPEYTWFKDGI VMPTNPKSTRAFSNSS YVLNPTTGELVFDPLASDT GEYSCEARNGYGTPMSTNAV RMEAVERNVGVIVAA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

