

Human Frizzled-1/FZD1 protein (Recombinant) (C-hFc&His) (STJP001070)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8321
Gene Symbol	FZD1
UniProt ID	FZD1_HUMAN
Immunogen Region	Gln73-His253
Immunogen Sequence	QAAGQGPGQGPGGQQPPPP PQQQQSGQQYNGERGISVPD HGYCQPISIP LCTDIAYNQT IMPNLLGHTNQEDAGLEVHQ FYPLVKVQCSAELKFFLC SM YAPVCTVLEQALPPCRSLCE RARQGCEALMNKFGFQWPD T LKCEKFPVHGAGELCVGQNT SDKGTPTPSLLPEFWTSNPQ H

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

