

Human HCVADR/CXADR protein (Recombinant) (C-His) (STJP001071)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1525
Gene Symbol	CXADR
UniProt ID	CXAR_HUMAN
Immunogen Region	Leu20-Gly237
Immunogen Sequence	LSITTPPEEMIEKAKGETAYL PCKFTLSPEDQGPLDIEWLI SPADNQKVDQVILYSGDKI YDDYYPD LKGRVHFTSNDLK SGDASINVTNLQLSDIGTYQ CKVKKAPGVANKKIHLLVVLV KPSGARCYYVDGSEEIGSDFK IKCEPKEGSLPLQYEWQKLS DSQKMPTSWLAEMTSSVISV KNASSEYSGTYSCTVRNRVG SDQCLLRLNVVPPSNKAG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

