

Human VCAM-1/CD106 protein (Recombinant) (C-His) (STJP001072)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7412
Gene Symbol	VCAM1
UniProt ID	VCAM1_HUMAN
Immunogen Region	Phe25-Glu698
Immunogen Sequence	FKIETTPESTRYLAQIGDSVS LTCSTTGCESPFFSWRTQID SPLNGKVTNEGTTSTLTMNP VSGNEHSYLCTATCESRKL EKGIQVEIYSFPKDPEIHLG GPLEAGKPITVKCSVADVYP FDRLEIDLKGDHLMKSQEF LEDADRKSLETKSLEVTFTP VIEDIGKVLVCRAKLHIDEM DSVPTVRQAVKELQVYISPK NTVISVNPSTKLQEGGVSMTM TCSSEGLPAPEIFWSSKLLD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

