

Human CD47 protein (Active Recombinant) (C-hFc&His) (STJP001073)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
UniProt ID	CD47_HUMAN
Immunogen Region	Gln19-Pro139
Immunogen Sequence	QLLFNKTKSVEFTFCNDTVV IPCFVTNMEAQNTTEVYVKW KFKGRDIYTFDGALNKSTVP TDFSSAKIEVSQLLKGDASL KMDKSDAVSHTGNYTCEVTE LTREGETIIEELKYRVVSWFS P

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

