

Human IFN-lambda 1/IL-29 protein (Recombinant) (C-His) (STJP001077)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	282618
Gene Symbol	IFNL1
UniProt ID	IFNL1_HUMAN
Immunogen Region	Gly20-Thr200
Immunogen Sequence	GPVPTSKPTTTGKGCHIGRF KSLSPQELASFKKARDALEE SLKLKNWSCSSPVFPGNWDL RLLQVRERPVALEAELALTL KYLEAAAGPALEDVLDQPLH TLHHILSQLQACIQPQPTAG PRPRGRLHHWLHRLQEAPKK ESAGCLEASVTFNLFRLTR DLKYVADGNLCLRTSTHPES T

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

