

Human PVR/CD155 protein (Recombinant) (C-hFc&His) (STJP001080)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5817
Gene Symbol	PVR
UniProt ID	PVR_HUMAN
Immunogen Region	Asp28-Asn343
Immunogen Sequence	DVVVQAPTQVPGFLGDSVTL PCYLQVPNMEVTHVSQLTWA RHGESGSMVAFHQQTGQPSYS ESKRLEFVAARLGAELRNAS LRMFGLRVEDEGNYTCLFVT FPQGSRSVDIWLRLAKPQN TAEVQKVQLTGEPVPMARCV STGGRPPAQITWHSDLGGMP NTSQVPGFLSGTVTITSLWI LVPSSQVDGKNVTCKVEHES FEKPQLLTVNLTVYYPPEVS ISGYDNNWYLGQNEATLTC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

